

# EMBRYO TRANSFER SET

## Set For Difficult Embryo Transfer

### Product description

- Guided catheter (polypropylene), with a sliding ring with centimetric graduations.
- Malleable removable mandrel.
- Transfer catheter (polyurethane) with centimetric graduations on the proximal part.
- The proximal (luer) end of the guiding catheter has been specifically designed for making the transfer catheter insertion easy.

### When to use?

- Embryo transfer in uterus, especially in the case of difficult transfers.
- Intrauterine insemination.

### How to use?

Before use, wash the embryo transfer set carefully with an appropriate sterile culture medium.

- With the help of centimetric graduations, position the sliding ring (guiding catheter) according to the depth of the uterus (Uterus depth should be measured prior to procedure).
- Give the guiding catheter-mandrel a shape, if necessary.
- Cross the cervical canal up to the internal orifice.
- Withdraw the mandrel.
- Insert the transfer catheter into the guiding catheter and push it up until the first centimetric graduation reaches the guiding catheter entrance. Distal ends of both catheters should match properly.
- Then push the transfer catheter with the help of centimetric graduations and proceed to transfer the embryo.
- Very slowly, withdraw the two catheters simultaneously.



Product Code	Product	Length		I.D.	O.D.	Sterilisation
EBT00118-3P	Embryo Transfer Set	18 cm	Guiding Catheter	1.7	2.3	γ irradiation
			Transfer Catheter	0.7	1.5	
EBT00123-3P	Embryo Transfer Set	23 cm	Guiding Catheter	1.7	2.3	γ irradiation
			Transfer Catheter	0.7	1.5	

Transfer catheter exceeds the guiding catheter by over 5 cm

Marketed in India by:

**BONNE SANTÉ**  
Translating Technology. Transforming Lives

Manufactured by: Laboratoires Ellios Bio Tek, PARIS, FRANCE.

