

# ELLIOCATH

## Curved Transfer Catheter

(Metal Free)

### Product Description

- Curved transfer catheter makes passage through the cervical canal easier to manage.
- Arrows located on the grip of ELLIOCATH® indicate the direction of the curve of the catheter.

### When to Use?

- Embryo transfer.

### How to Use?

Before use, wash the catheter carefully with an appropriate sterile culture medium.

### Embryo Transfer

- Place the embryo in ELLIOCATH® with a small volume of culture medium.
- Position the sliding ring according to the uterus depth. (Depth should be measured prior to the procedure).
- Advance ELLIOCATH® to a distance of 1 cm from the uterine fundus.
- To facilitate cervical canal passage, push the catheter with a gentle rotational movement.
- Transfer embryo.
- Slowly withdraw the catheter.



Product Code	Product	I.D.	Length	Sterilisation
NS1280545	ELLIOCATH®	1.1 mm	20.5 cm	γ irradiation

Marketed in India by:

**BONNE SANTÉ**  
Translating Technology. Transforming Lives

Manufactured by: Laboratoires Ellios Bio Tek, PARIS, FRANCE.

