

OOCYTES ASPIRATION NEEDLE

Back Cut Bevel

Product description

- Single lumen needle.
- Back-cut bevel (17G / 1.4mm).
- Non-coring bevel.
- Two laser echo markings.
- Plastic stylet.
- Large ergonomic handle.
- Tactile mark indicating bevel position.
- 4 tabs for traceability.

When to use?

- Echo guided oocytes pick up.
- Ovarian cyst aspiration.

How to use?

Before use, wash the needle carefully with an appropriate sterile culture medium.

- Adapt a syringe to the luer tip.
- Proceed according to the selected route.



Product Code	Product	O.D.	Length	Sterilisation
EBT00030-BC	Ovum aspiration needle with stylet	17 G (1.4 mm)	30 cm (active length)	γ irradiation
EBT00030-BC-O	Ovum aspiration needle without stylet	17 G (1.4 mm)	30 cm (active length)	γ irradiation

Marketed in India by:

BONNE SANTÉ
Translating Technology. Transforming Lives

Manufactured by: Laboratoires Ellios Bio Tek, PARIS, FRANCE.

